

The following compounds show no cross-reactivity when tested with **AccuSign® BZO** at a concentration of 100 µg/mL. (Table 3.)

**Table 3. Non Cross-Reacting Compounds**

Acetaldehyde	Ecgonine methyl-ester	Methaqualone
4-Acetamidophenol	(+) Ephedrine	Methoxyphenamine
Acetaminophen	(±) Ephedrine	(±) 3,4-Methylene-dioxyamphet-amine
Acetone	(-) Ephedrine	
Acetophenetidin (Phenacetin)	(-) Ψ Ephedrine	(±) 3,4-Methylene-dioxymetham-phetamine
N-Acetylprocainamide	Epinephrine	Methylphenidate
Acetylsalicylic acid	Erythromycin	Methyprylon
Albumin	β-Estradiol	Morphine
Aminopyrine	Estrone-3-sulfate	Morphine-3-β-D-glucuronide
Amitypytline	Ethyl-p-amino-benzoate	Nalidixic acid
Amobarbital	Fenopropfen	Nalorphine
Amoxapine	Furoxide	Naloxone
Amoxicillin	Gentisic acid	Naltrexone
D,L-Amphetamine	Guaiacol glycerol ether	Naproxen
L-Amphetamine		Niacinamide
Ampicillin	Glucose	Nifedipine
Apomorphine	Glucuronide	Norcodein
Ascorbic acid	Glutethimide	Norethindrone
Aspartame	Guaifenesin	Noroxymorphone
Atropine	Hemoglobin	D-Norpropoxyphene
Benzilic acid	Hippuric acid	(-) Norpseudo-ephedrine
Benzocaine	Hyalalazine	Noscapine
Benzoic acid	Hydrochlorothiazide	Nylidrin
Benzoylcegonine	Hydrocodone	D,L-Octopamine
Benzphetamine	Hydrocortisone	Oxalic acid
Bilirubin	Hydromorphone	Oxolinic acid
Butabarbital	O-Hydroxyhippuric acid	Oxycodone
Cannabidiol	3-Hydroxytyramine	Oxymetazoline
Cannabinol	Ibuprofen	Oxymorphone
Chloralhydrate	Imipramine	Papaverine
Chloramphenicol	Iproniazid	Penicillin-G
Chlorothiazide	(-) Isoproterenol	Pentazocaine
Chlorpheniramine	Isoxsuprine	Pentobarbital
Chlorpromazine	Ketamine	Perphenazine
Chloroquine	Ketones	Phencyclidine
Cholesterol	Cocaine hydrochloride	Phendimetrazine
Clomipramine	Ketoprofen	Phenelzine
Clonidine	Labetalol	β-Phenethylamine
Cocaine hydrochloride	Levorphanol	Phenobarbital
Codeine	Lidocaine	Phenothiazine
Cortisone	Loperamide	Phentermine
(-) Cotinine	Loxapine succinate	Phenoin
Creatinine	Lysergic acid diethylamide	Phenylbutazone
Deoxycorticosterone	Maprotiline	L-Phenylephrine
Dextromethorphan	Melanin	D,L-Phenylpropanol-amine
Dextropropoxyphene	Meperidine	Prednisolone
Diclofenac	Meprobamate	Prednisone
Diethylpropion	Methadone	Procaine HCl
Diflunisal	D-Methamphet-amine	Promazine
Digoxin	Doxylamine	Promethazine
Diphenhydramine	p-Hydroxymeth-amphetamine	
Domperidone		
Doxylamine		
Ecgonine hydrochloride		

D,L-Propranolol  
Propiomazine  
D-Propoxyphene  
D-Pseudoephedrine  
L-Pseudoephedrine  
Quinidine  
Quinine  
Rantidine  
Salicylic acid  
Secobarbital  
Serotonin  
Sodium chloride  
Sulfamethazine


















Sulindac  
Tetracycline  
Tetrahydrocortisone  
Δ<sup>8</sup>-THC  
Δ<sup>9</sup>-THC  
11-nor-Δ<sup>9</sup>-THC-9-COOH  
Tetrahydrozoline  
Thebaine  
Thiamine  
Thioridazine  
D,L-Thyroxine  
Tolbutamide

Triamterene  
Trifluoperazine  
Trimethoprim  
Trimipramine  
Tryptamine  
D,L-Tryptophan  
Tyramine  
D,L-Tyrosine  
Uric acid  
Verapamil  
Zomepirac

## References

- Hawks RL, Chiang CN, eds. *Urine Testing for Drugs of Abuse*. Rockville, MD: National Institute on Drug Abuse (NIDA), Research Monograph 73;1986.
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*, 2nd Ed., Davis, CA: Biomedical Publ.;1982.
- Greenblatt DJ, Shader RI. *Benzodiazepines in Clinical Practice*. New York: Raven Press;1974.

## Symbols Key


	Manufactured by
	CE Mark
	Authorized Representative
	In Vitro Diagnostic Medical Device
	Catalog Number
	Consult Instructions for Use
	Batch Code
	“Use By” date in year-month-day format
	Temperature Limitation
	Contains sufficient for <n> tests
	Do not reuse
	Contents
	Test Device
	Transfer Pipette
	Instructions for Use
	One-step immunochromatographic Assay for the Detection of Drugs of Abuse in Urine
	Benzodiazepines Test

AccuSign® is a Registered Trademark of Princeton BioMeditech Corporation.

Patent No.: 5,559,041



© 2000 PBM  
Printed in U.S.A.  
Revised May 2012  
P-5846-F 0530BL

   
MTPromedConsultingGmbH  
Altenhofstrasse 80  
66386 St. Ingbert  
Germany  
+49-68 94-58 10 20

Manufactured by  
  
Princeton BioMeditech Corporation  
4242 US Hwy 1  
Monmouth Jct., NJ 08852, U.S.A.  
1-732-274-1000 www.pbmc.com

P-5846-F

# AccuSign® BZO

## One-Step Benzodiazepines Test

For *In Vitro* Use Only

### Simple One-Step Immunoassay for the Qualitative Detection of Benzodiazepines and/or their Metabolites in Human Urine

## PBM

Catalog No. DOA-207-35 35 Test Kit  
DOA-207-10 10 Test Kit

## Intended Use

**AccuSign® BZO** is a simple, one-step immunoassay intended for use in the qualitative detection of benzodiazepines in human urine with a cutoff at 300 ng/mL for oxazepam.

*The AccuSign® BZO test provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography, mass spectrometry (GC/MS) is the preferred confirmatory method. Other chemical confirmatory methods are available. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.*<sup>1</sup>

## Summary and Explanation

Benzodiazepines are a class of frequently prescribed central nervous system (CNS) depressants which include widely used drugs such as chlordiazepoxide, diazepam, and oxazepam. They have medically useful properties, including antianxiety, sedative, anti-convulsant, and hypnotic effects.<sup>2</sup> They are taken orally or sometimes by injection, and have a low potential for physical or psychological dependence. Benzodiazepines induce drowsiness and muscle relaxation. Their use can also result in intoxication, similar to drunken behavior except without evidence of alcohol use, and the loss of inhibitions. Chronic abuse can result in addiction and tardive dyskinesia (involuntary muscle movements of the face, limbs, and trunk). Overdose can result in coma and possible death. Withdrawal syndrome includes anxiety, insomnia, tremors, delirium, and convulsions.

The effects of benzodiazepine last 4–8 hours. The different benzodiazepines are absorbed at different rates, and the timing of their psychoactive effects varies with the absorption rate. The drugs are excreted in the urine primarily as the parent compounds or as oxazepam glucuronide, an inactive metabolite (in the case of chlordiazepoxide and diazepam) and are detectable for 1–2 days. Oxazepam is detectable in the urine for up to 7 days.<sup>2,3</sup>

## Principle

The **AccuSign® BZO** test uses solid-phase chromatographic membrane immunoassay technology for the qualitative detection of benzodiazepine. The test is based on the principle of the highly specific immunochemical reactions between antigens and antibodies which are used for the analysis of specific substances in biological fluids. The test relies on the competition between the drug conjugates and the drugs which may be present in the urine sample, for binding to antibodies. In the test procedure, a sample of urine is placed in the Sample well of the device and is allowed to migrate upward. If the drug is present in the urine sample, it competes with the drug conjugate bound to the dye, for the limited antibodies immobilized on the membrane. If the level of drug or drug metabolite is above the cutoff level, the drug will saturate the antibodies, thus inhibiting the binding of the dye coated with drug conjugates to the antibodies on the membrane. This prevents the formation of a line on the membrane. Therefore, a drug-positive urine sample will not generate a line at Test (T) position in the Result window, indicating a positive result from positive drug competition. A negative urine sample will generate a line at Test position in the Result window, indicating a negative result from an absence of competition with free drugs.

In addition to the line that may appear at the Test position in the Result window, a Control line must appear at the Control (C) validation position in the Result window to confirm the viability of the test. This Control line should always be seen if the test is conducted properly. This works as a procedural control, confirming that proper sample volume was used and the reagent system worked. If insufficient sample volume is used, there may not be a Control line, indicating that the test is invalid.


## Materials Provided

The **AccuSign® BZO** test kit contains all the reagents necessary to perform the tests.

- AccuSign® BZO** device. The test device contains a membrane strip coated with polyclonal anti-benzodiazepine antibody and a pad containing drug-dye conjugate in a protein matrix.
- Disposable specimen dispensers.
- Instructions for use.

## Precautions

- For *in vitro* diagnostic use only.
- Avoid cross contamination of urine samples by using a new urine specimen container and dropper for each urine sample.
- The test kit does not contain any HIV or hepatitis infective components.
- Urine specimens are potentially infectious. Proper handling and disposal methods should be established according to good laboratory practices.
- The **AccuSign®** device should remain in its original sealed pouch until ready for use. Do not use the test if the pouch is damaged or the seal is broken.
- Do not use the test kit after the expiration date.



### Storage and Stability

The **AccuSign® BZO** test kit should be stored at 2–30°C (35–86°F) in the original sealed pouch. The expiration dating was established under these storage conditions.

### Specimen Collection and Preparation

Approximately 110 µL of urine sample is required for each test. Fresh urine specimens do not require any special handling or pretreatment. Specimens should be collected in a clean glass or plastic container. If testing will not be performed immediately, specimens should be refrigerated (2–8°C) or frozen. Specimens should be brought to room temperature before testing.

Specimens containing a large amount of particulate matter may give inconsistent test results. Such specimens should be clarified by centrifuging or allowing to settle before testing.

### Test Procedure

The test procedure consists of adding the urine sample to the Sample well of the device and watching for the appearance of colored lines in the result window.

#### Test Protocol

1. For each test, open one **AccuSign® BZO** pouch and label the **AccuSign®** device with the patient ID.
2. Holding the dropper vertically, dispense 3 drops (110 µL) of the urine sample into the Sample well (S).
3. Read the result after 3 minutes, but within 10 minutes of sample addition.

### Interpretation of Results

**Negative:** The appearance of a reddish-purple Control line (C) and a line next to T indicates a negative test result; i.e., no drug above the cutoff level has been detected. The color intensities of the Control line and the Test line may not be equal. *Any faint line in the Result window, visible in 10 minutes, should be interpreted as negative. A negative test result does not indicate the absence of drug in the sample; it only indicates the sample does not contain drug above the cutoff level in qualitative terms.*

**Positive:** The appearance of only a reddish-purple Control line and no distinct line next to T indicates the test result is positive for BZO (i.e., the specimen contains BZO at a concentration above the cutoff level). *A positive test result does not provide any indication of the level of intoxication or urinary concentration of the drug in the sample; it only indicates the sample contains drug above the cutoff level in qualitative terms.*

**Invalid:** A distinct Control line (C) should always appear. The test is invalid if no Control line forms at the C position. Such tests should be repeated with a new **AccuSign® BZO** test device.

### Limitations

- The test is designed for use with unadulterated urine only.
- There is a possibility that factors such as technical or procedural errors, as well as other substances in the urine sample which are not listed in Table 3 below, may interfere with the test and cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results. If adulteration is suspected, the test should be repeated with a new sample.
- The test result read after 10 minutes may not be consistent with the original reading obtained within the 10 minute reading period. The test must be read within 10 minutes of sample application.

### User Quality Control

**Internal Control:** Each **AccuSign®** test device has built-in controls. The Control line is an internal positive procedural control. A distinct reddish-purple Control line should always appear at the C position, if the test procedure is performed properly, an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the test reagents at the control line and the conjugate-color indicator are reactive. In addition, if the test has been performed correctly and the device is working properly, the background in the result window will become clear and provide a distinct result. This may be considered an internal negative procedural control.

The positive and negative procedural controls contained in each **AccuSign®** test device satisfy the requirements of testing a positive control and a negative control on a daily basis. If the Control line does not appear at the Control position, the test is invalid and a new test should be performed. If the problem persists, contact PBM for technical assistance.

**External Control:** External controls may also be used to assure that the reagents are working properly and that the assay procedure is followed correctly. It is recommended that a control be tested at regular intervals as good laboratory testing process. For information on how to obtain controls, contact PBM's Technical Services.

### Expected Values

**AccuSign® BZO** is a qualitative assay. The amount of benzodiazepines and/or their metabolites in the urine cannot be estimated by the assay. The assay results distinguish positive from negative samples. Positive results indicate the samples contain benzodiazepine metabolites above the cutoff concentration.

### Performance Characteristics

The **AccuSign® BZO** test has been shown to detect oxazepam at an average cutoff of 300 ng/mL in urine. The test also detects other benzodiazepines listed below at the minimum concentrations indicated (Table 2).

The accuracy of **AccuSign® BZO** was evaluated in comparison to a commercially available immunoassay (Syva® EMIT® II) at a cutoff of 300 ng/mL. A total of 223 samples was tested by both procedures. The overall accuracy of the test was 98.7%, as shown below. (Table 1.)

**Table 1. Accuracy: Comparison of AccuSign® BZO with Syva® EMIT® II**

		Syva® EMIT® II (BZO)		
		Positive	Negative	TOTAL
<b>AccuSign® BZO</b>	Positive	84	2	86
	Negative	1	136	137
<b>TOTAL</b>		85	138	223
		Relative Sensitivity		Relative Specificity
<b>AccuSign® BZO</b>		98.8% (84/85)		98.6% (136/138)

Discrepant samples for BZO were analyzed by GC/MS. The one false-negative sample contained the drug at a level of 344 ng/mL, while the two false-positive samples showed 274 and 289 ng/mL.

In a separate study, **AccuSign® BZO** was evaluated against 27 specimens confirmed as positive by GC/MS. The range of drug values was 312 to 610 ng/mL. The results demonstrate the excellent correlation of **AccuSign® BZO** with GC/MS.

### Precision and Accuracy

The precision of the **AccuSign® BZO** assay was determined by carrying out the test with serially diluted standard drug solutions using 3 lots of products on 3 different dates. Ninety-eight percent (98%) of the spiked samples containing oxazepam concentrations 25% over the cutoff level (i.e., 375 ng/mL) showed positive results.

### Distribution of Random Error:

Twenty (20) blind samples prepared by spiking various concentrations of drug were separately tested by two operators. The test results from the two operators showed complete agreement.

### Reproducibility

The reproducibility of the test results of the **AccuSign® BZO** assay was examined at three different sites using a total of 15 blind controls, consisting of 5 negative samples, 5 moderately positive samples (600 ng/mL oxazepam), and 5 strongly positive samples (1,200 ng/mL oxazepam). The results obtained at these three sites with these controls demonstrated 100% agreement with each other.

### Specificity

Compounds that are detected by the **AccuSign® BZO** test are listed below (Table 2). The specificity of **AccuSign® BZO** was determined by adding various drugs and drug metabolites to drug-negative urine specimens and testing with the **AccuSign® BZO** test kit. The results are expressed in terms of the concentration required to produce a positive result.

**Table 2. Specificity**

Compound	Concentration (ng/mL)
Alprazolam	100,000
Bromazepam	1,250
Chlordiazepoxide	500
Clobazam	>100,000
Clonazepam	30,000
Clorazepate dipotassium	2000
Delorazepam	1,500
N-Desalkylflurazepam	2,500
Diazepam	10,000
Estazolam	>100,000
Flunitrazepam	>100,000
7-amino-flunitrazepam	1,500
a-Hydroxyalprazolam	100,000
a-Hydroxytriazolam	10,000
Lorazepam	2,500
Lormetazepam	25,000
Medazepam	10,000
Midazolam	25,000
Nitrazepam	100,000
Nordiazepam(N-Desmethyldiazepam)	7,500
Oxazepam	300
Prazepam	>100,000
Temazepam	6,000
Triazolam	>100,000