

Canine Heartworm Antigen Test Kit

BioSign® CHW

Rapid test for the detection of adult *Dirofilaria immitis* antigen in canine serum, plasma or whole blood.



- Results in 10 Minutes
- Whole blood, serum or plasma sample
- Room temperature storage
- Visual read or optional instrument read with DXpress Reader
- Sensitivity 97.9% and Specificity 97.5%

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
Canine Heartworm Antigen Test Kit BioSign CHW	30	BSP-604
Canine Heartworm Antigen Test Kit BioSign CHW	10	BSP-604-10

 Proudly made in the USA

For more information or to place an order, contact:



A PBM Group Company
 85 Orchard Road • Skillman • New Jersey • 08558
 800.526.2125 • 732.246.3366 • Fax: 732.246.0570