

Cultura[®] M

Cultura[®] M is a small, durable and reliable incubator suitable for the incubation of dip paddles.



Cultura[®] M has been used for 10 years in the fields of medical diagnostics, food bacteriology, hygiene monitoring, and industrial process testing.

SPECIFICATIONS

Temperature range: 25 - 45°C	Case: Not easily flammable PC & ABS		
Temperature stability: +/- 1°C	Door: PC transparent		
Voltage: 12 V / 110 V / 230 V	Thermometer: up to 60°C, in special tube		
Frequency: 50 Hz	Base cover: ABS		
Power output: 26 W / 0.2 A	Warranty: 2 years from initial installation		
Protective Class: II	Weight: 2.4 lbs (1.1 kg)		
Dimensions:	W	H	D
External	12.2 in (310 mm)	6.1 in (155 mm)	6.6 in (168 mm)
Internal	8.7 in (220 mm)	4.7 in (120 mm)	5.9 in (150 mm)

For more information about
Cultura[®] M Incubator
or to place an order, contact:

A PBM Group Company

71 Veronica Avenue
Somerset, New Jersey 08873
800.526.2125 732.246.3366
Fax: 732.246.0570
www.lifesignmed.com

Cultura[®] M is a registered trademark of Almedica AG.
LifeSign is a registered trademark of Princeton BioMeditech Corporation.

This CE Mark indicates that this equipment complies with the following directive:

- EMC 89/336/EEC
- DIN EN 61010-4
- 61010-2010



ORDERING INFORMATION

PRODUCT	CATALOG #
Cultura [®] M Incubator	05866